

Nr.crt	Descriere	Pret unitar Euro
1.	Microscop endodontic OPMI pico <ul style="list-style-type: none">- marire in 5 trepte- lentila obiectiv $f=250\text{mm}$ cu inel focusare- tub binocular inclinat- stand podea S100- sursa iluminare Ha cu 2 becuri 12V/100W- schimbare rapida a becului ars- fibra optica- maner manevrare	10.000,00



Accesorii microscop:

2. Sistem video integrat, PAL, 1 CCD, monitor plat TFT (camera video analogica) 4.700,00



3. Sistem video HD integrat 720 p (fara monitor) 5.000,00



4. Sistem video HD integrat 1080 p (fara monitor) 6.400,00



5. Separator optic si adaptor FlexioStill
Pentru camera foto/video digitala (fara camera !!) 3.950,00



6. Separator optic si adaptor camera foto digitala SLR
(fara camera si adaptor T2!!)

3.850,00



7. Sistem optic MORA, cu miscare laterala si inainte-inaoi
(cost suplimentar !!)

5.950,00



© Carl Zeiss

8. Tub binocular cu miscare inainte-inaoi 180 grade
(cost suplimentar !!)

2.400,00



9. Sursa iluminare Xenon cu rezerva
(cost suplimentar !!)

3.900,00



Preturile nu contin TVA.

Modalitate de plata: in Lei, la curs BNR, 40% in avans, 60% in 5 zile de la livrare.

Termen de livrare: 6 saptamani de la plata avansului.

Se asigura: livrare, instalare, instruire si 12 luni garantie de la instalare.

- | | |
|---|-----------|
| 1. Microscop operator dental OPMI PROergo/S7 | 34.600,00 |
| <ul style="list-style-type: none"> - optica corectie apocromata - sistem Varioskop (modificare distanta de lucru 200 ó 415 mm) motorizat - zoom motorizat - filtre optice verde si portocaliu - tub binocular inclinabil - manere multifunctionale - stand podea S7 - frane electromagnetice pt. pozitionare microscop - posibilitate sustinere 14 kg - sistem balansare - iluminare fibra optica - sursa lumina Ha - 2 x bec Ha cu schimbare automata | |



- | | |
|--|----------|
| 2. Sistem video integrat, PAL, 1 CCD si monitor plat
(camera video analogica) | 5.900,00 |
| 3. Sistem optic cu miscare laterala, separator optic si adaptor
camera foto digitala (fara camera !!) | 5.850,00 |
| 4. Sistem optic cu miscare laterala | 2.900,00 |
| 5. Sursa iluminare Xenon cu rezerva
(cost suplimentar !!) | 3.950,00 |

Preturile nu contin TVA.

Modalitate de plata: in Lei, la curs BNR, 40% in avans, 60% in 15 zile de la livrare,
Termen de livrare: 6 saptamani de la plata avansului.

Se asigura: livrare, instalare, instruire si 12 luni garantie de la instalare.

Nr.crt	Descriere	Pret unitar Euro
1.	Prisme chirurgicale cu banda de cap EyeMag Pro S	1.950,00



2.	Prisme chirurgicale cu ochelari cu rama de titan EyeMag Pro F	1.700,00
----	---	----------



3.	Sistem iluminare cu halogen pt. prisme chirurgicale	1.800,00
----	---	----------

4.	Lupe EyeMag Smart	1.280,00
----	-------------------	----------



5. Sistem de iluminare cu led

990,00



Preturile nu contin TVA.

Modalitate de plata: in Lei, la curs BNR, 40% in avans, 60% in 15 zile de la livrare,

Termen de livrare: 3 saptamani de la plata avansului.

Se asigura: livrare, instruire si 12 luni garantie de la livrare.