

# X-SMART™ EASY

THE CORDLESS ENDO MOTOR  
FOR THE GENERAL PRACTITIONER



**DENTSPLY**  
MAILLEFER

# X-SMART™ EASY

CORDLESS – LIGHTWEIGHT – SIMPLE TO USE

## IMPROVED FIELD OF VISION & ACCESS

- Slim contra-angle head, adjustable in 6 positions

## FREEDOM OF MOVEMENT

- Cordless
  - 1 hour continuous use without recharge
  - Full recharge time in less than 2 hours
  - Automatic power-off after 3 minutes of non-use
- Lightweight

## FLEXIBILITY

Complete set of torque and speed settings

## SAFETY

Auto reverse at torque limit



Forward button

Reverse button

Press any button to turn the unit on

LCD Screen: display of the selected parameters

### Keyboard for selection of settings:

- Speed: 100 to 800 rpm
- Torque: 0,6 to 4 Ncm
- Auto reverse: On/Off

**X-SMART™ EASY**  
Motor + contra-angle  
REF A 1022

**Contra-angle 16:1**  
REF A 1006

**X-SMART™ Spray**  
Lubrication spray, 500ml  
REF A 1007-8

**X-SMART™ Cleaning Adaptor**  
Contra-angle adaptor for automatic  
cleaning and maintenance units  
(F-type to Midwest 4-hole)  
REF A 1007-6

