

thermaprep[®] 2

3-dimensional heat
for a 3-dimensional obturation



- 3-dimensional heat
- Indicator lights and audible alerts
- Slow release obturator holders
- Cleaning mode

thermaprep® 2

3-dimensional heat for a 3-dimensional obturation

3-dimensional heat

Obturator holders are heated in seconds in single and dual heating modes. 3-Dimensional heat is evenly maintained in the heating chamber until the obturator is ready to use, making it ideal for multi-canal procedures.

Indicator lights and audible alerts

The oven uses both LED indicator lights and audible alerts to indicate status changes.



Slow release obturator holders

Obturator holders are designed to allow stable vertical entry and exit of the heating chamber.

Cleaning mode

The cleaning mode raises the temperature of the oven and softens gutta-percha inside the heating chamber for removal.

- The oven heats up obturators in 20 to 49 seconds, depending on the size of the obturator
- A built-in timer audibly signals when the filling material reaches the right temperature
- The obturator is kept at the right temperature for 90 seconds.

thermaprep® 2 oven

Contents:

- Thermaprep® 2 oven
 - Cleaning brush
 - Additional faceplate stickers*
- EU 230V REF A 0176 230 002 00
 UK 230V REF A 0176 230 002 33
 115V REF A 0176 115 002 00

thermaprep® 2 oven kit for ProTaper® Universal

Contents:

- Thermaprep® 2 oven
 - ProTaper® universal obturators F1 x6 REF A 1411
 - ProTaper® universal obturators F2 x20 REF A 1411
 - ProTaper® universal obturators F3 x6 REF A 1411
 - Cleaning brush
 - Additional faceplate stickers*
- EU 230V REF A 176C 230 002 00
 UK 230V REF A 176C 230 002 33
 115V REF A 176C 115 002 00

thermaprep® 2 oven kit for Thermafil®

Contents:


- Thermaprep® 2 oven
 - Thermafil® verifiers assortment 20-45 REF A 0175
 - Thermafil® obturators posterior kit x20 REF A 0167
 - Cleaning brush
 - Additional faceplate stickers*
- EU 230V REF A 176A 230 002 00
 UK 230V REF A 176A 230 002 33
 115V REF A 176A 115 002 00

thermaprep® 2 oven kit for wave•one™

Contents:

- Thermaprep® 2 oven
 - wave•one™ obturators small x6 REF A 1700
 - wave•one™ obturators primary x20 REF A 1700
 - wave•one™ obturators large x6 REF A 1700
 - Cleaning brush
 - Additional faceplate stickers*
- EU 230V REF A 176B 230 002 00
 UK 230V REF A 176B 230 002 33
 115V REF A 176B 115 002 00

* In addition to the Thermafil® obturators, the Thermaprep®2 oven works with ProTaper® Universal obturators, wave•one™ obturators and GT® SeriesX obturators. The additional faceplate stickers are included in all versions in order to allow the adaptation of the oven to the chosen system.

Distributed by:  Dentsply Tulsa Dental Specialties
 608 Rolling Hills Drive
 Johnson City, TN 37604

Made in the Netherlands