



Nr.crt	Descriere	Pret unitar Euro
1.	Prisme chirurgicale cu banda de cap <b>EyeMag Pro S</b>	2.300,00



2.	Prisme chirurgicale cu ochelari cu rama de titan <b>EyeMag Pro F</b>	2.100,00
----	----------------------------------------------------------------------	----------



3. **Lupe EyeMag Smart**

1.500,00



4. **Sistem de iluminare LED EyeMag Light II**

1.400,00



Preturile nu contin TVA.

Termen de livrare: 6 saptamani de la plata avansului.

Se asigura: livrare, instalare, instruire si 12 luni garantie de la instalare.

The ZEISS logo is displayed in white, bold, sans-serif capital letters on a blue rectangular background with a curved bottom edge.

Nr.crt	Descriere	Pret unitar Euro
1.	<b>Microscop endodontic OPMI pico</b> <ul style="list-style-type: none"><li>- marire in 5 trepte</li><li>- lentila obiectiv <math>f=250\text{mm}</math> cu inel focusare</li><li>- tub binocular inclinat</li><li>- stand podea S100</li><li>- sursa iluminare Ha cu 2 becuri 12V/100W</li><li>- schimbare rapida a becului ars</li><li>- fibra optica</li><li>- maner manevrare</li></ul>	12.500,00



#### Accesorii microscop:

2.	Sistem video HD integrat 720 p (fara monitor)	5.500,00
----	-----------------------------------------------	----------




3. Sistem video HD integrat 1080 p (fara monitor) 6.900,00



4. Separator optic si adaptor FlexioStill  
Pentru camera foto/video digitala (fara camera !!) 4.900,00



5. Separator optic si adaptor camera foto digitala SLR  
(fara camera si adaptor T2!!) 3.900,00



6. Sistem optic MORA, cu miscare laterala  
necesita cuplare cu tub binocular inclinabil sau pliabil  
(cost suplimentar !!) 4.900,00




7. Tub binocular inclinabil cu miscare inainte-inapoi 180 grade (cost suplimentar !!) 2.900,00



8. Tub binocular pliabil, cu camp vizual largit, sistem aditional de marire cu factor 1,5 x integrat in tubul binocular, inel rotatie tub binocular integrat cu piedica (cost suplimentar !!) 5.900,00



9. Sursa iluminare Xenon cu rezerva (cost suplimentar !!) 4.500,00



10. Sursa iluminare LED (cost suplimentar !!) 1.900,00



11. **Sistem varioscopic** 3.300,00  
permite modificarea continua a distantei de lucru 200-300 mm  
(cost suplimentar !!)



12. **Set lentile pentru protectie lentila obiectiv** 400,00



Preturile nu contin TVA.

Termen de livrare: 6 saptamani de la plata avansului.

Se asigura: livrare, instalare, instruire si 12 luni garantie de la instalare.